



Gel Warmer

GW-210/GW-220

The GW-200 gel warmer can be used to:

- Increase gel liquidity, easy to use
- Adjust gel temperature, fell comfortable
- Improve gel use quality effect

The GW-200 gel warmer can be used when using ultrasound machine, ECG, fetal doppler, fetal monitor, physiatrics instrument and lithotripter to human and animal.



CE



Features

- Displaying parameters with LED
- Adjustable temperature range 34-40°C (for GW-220)
- Programmable timer to switch on and off (for GW-220)
- 360 degree surround heater unit
- Digital temperature control technology, accuracy: $\pm 1^{\circ}\text{C}$
- High energy efficiency
- Ergonomic 2 bottles design
- Built-in wall mountable construction