

LANCING DEVICES


Adjustable Automatic Lancing Devices



Assist diabetic patients to monitor their glucose levels, easily and safely with an adjustable dial, lancing depth can be set to an appropriate level for the individual. Any standard lancet can be used with these devices.



LANCING DEVICE



Patent




HOW TO USE **LANCING DEVICES**


Please read instructions carefully before use.
Keep hands clear at all time during use!



1 Quick, easy lancet insertion and disposal.



2 It is simple to remove and twist off lancet protective cap.



3 Five depth settings for optimum skin penetration.



4 Sliding barrel signals readiness with a soft click.



5 One-handed lancing with comfortable trigger button.

Reorder#	Type	Color	Depth	Print	Package
ZWCB-1105	LANCING DEVICE Normal	White gray / Blue	9	5	1/bag 200/case
ZWCB-1109	LANCING DEVICE Normal	White gray / Blue		9	