

HEALTH MONITORING



Saliva Alcohol Test Strip

One Step Saliva Alcohol Screening Test is a simple, accurate and convenient chemical assay for the rapid and qualitative detection of human blood alcohol levels by comparing the color developed at the end of the reaction period with the color chart provided.

Features

1. Fast and easy to use, no other reagent is needed
2. Semi-Quantitative results are read visually without the need for instrumentation
3. Sensitivity: 0.02% alcohol concentration in saliva.

Micro-Albumin Urine Rapid Test

When there is an abnormally high permeability for albumin in the renal glomerulus, albumin will leak into the urine. Usually, when the level is above 100 mg/L, it is considered serious micro-albuminuria; when it is over 300 mg/L, it is then called macro-albuminuria or albuminuria. The level of micro albumin can be used in many situations in clinical diagnosis realm, among them including: as an indicator of subclinical cardiovascular disease; as a marker of vascular endothelial dysfunction; as an important prognostic marker for kidney disease, such as diabetes mellitus and hypertension; and, as a risk factor for venous thromboembolism.

The Micro Albumin urine test kits are a chromatographic immunoassay for the semi-quantitative detection of micro-albumin in human urine specimen based on competitive inhibition principle. The test result can be used in many clinical diagnostic situations, especially as a prognostic marker of kidney damages.

Features

1. Fast and easy to use, no other reagent is needed
2. Semi-Quantitative results are read visually without the need for instrumentation
3. Positive Results: 0-20 mg/L; 20-50 mg/L; over 50 mg/L.

Ordering information

CAT. NO.	Description	Method	Packing	Specimen
ALC 711	Saliva Alcohol Rapid Test Strip	Chemical	25 Tests/Box; 160 Boxes/Carton	Saliva
ALC712	Saliva Alcohol Rapid Test Strip	GICA	25 Tests/Box; 160 Boxes/Carton	Urine
HSA 151	Micro-Albumin Urine Test Strip	GICA	100 Tests/Box; 48 Boxes/Carton	Urine
HSA 152	Micro-Albumin Urine Cassette	GICA	25 Tests/Box; 40 Boxes/Carton	Urine

* GICA: Lateral Flow Gold Immunochromatographic Assay; Chemical: Chemical colorimetric Assay;

** Carton Outer Size: 52*45*36CM. Gross weight about 14-15.5 KG/Carton.