

PD-668



PD Meter

Measurement Method	Cornea reflection light coincidence method
Range of data display	Binocular pupillary distance: 44-83mm
	Left or right pupillary distance: 22-41.5mm
Effective range of measurement	Binocular pupillary distance: 45-82mm
	Left or right pupillary distance: 22.5-41mm
Indication error	≤0.5mm
Rounding error	≤0.5mm
The scope of the focus point	30cm-∞
Power source	Voltage: DC 3V, Specification: AA battery, Quantity: 2 or 4 pieces
Time for automatic shut-off	1 minute after stopping operation
Dimensions	214mm(L) x 150mm(W) x 51mm(H)
Weight	0.6kg

BL-16

Motorized Table

Table size	57*38(CM)
Table elevation range	63 – 82 (CM)
Power	220V/50HZ or 110V/60HZ
Weight	20KGS
Max loading	70KGS

